Universal-Cleaner for the ultrasonic bath

Concentrate - mildly alkaline - with corrosion protection

EXAM-expertise No.: 5734/06

For cleaning and degreasing of machines and their components, tools, equipment and parts in industry, handicraft and laboratories. For cleaning of heavily soiled respirators (face masks) and protective work equipment.

TICKOPUR R 33 has been especially designed for cleaning in SONOREX ultrasonic cleaning units, as well as for cleaning by immersion, wiping or high pressure techniques. If used properly, a fast and gentle cleaning is achieved.

TICKOPUR R 33 removes:

common soiling, grinding-, polishing- and metals, also light metals, glass, ceramics, lapping-residues, oily and greasy residues, plastics, rubber, windows, glasses, E-filters, from dust, soot, ink etc. respirators (face masks) etc. **Properties** Material compatibility TICKOPUR R 33 has a high material compatibility. liquid Parts made of tin and zinc may be damaged. We highly concentrated recommend testing the resistance of sensitive for ultrasonic application materials before cleaning. • pH 9.9 at 1% The material compatibility of TICKOPUR R 33 in the • with corrosion protection ultrasonic bath on face masks of MSA Auer has been EXAM-expertise for ultrasonic application certified by an expertise of EXAM "centre of biodegradable competence" for breathing protection. Application Chemical and physical data Application with ultrasound appearance: clear, colourless liquid dosage: 3 - 5 % density: 1.13 g/ml sonication time: 1 - 10 minutes viscosity: 5.4 mPa·s 20 - 80 °C temperature: -6 to 100 °C temperature stability: Application without ultrasound shelf life: >6 years dosage: 1 - 10 % pH-value: 9.9 at 1 % immersion time: up to 12 hours Ecology 20 - 80 °C temperature: The contained surfactants are biodegradable Information for application according to Regulation (EC) No. 648/2004 on Make sure the surface of objects is completely detergents. Used solutions can be let into the immersed. Hollow objects have to be immersed in a sewerage after neutralisation to a pH-value of 6.5 to way that air can completely escape from the cavities. 9.5 and dilution with the quadruple amount of water if Concentration and time of application depend on the contaminations allow this. The national and local

type and degree of contamination. After cleaning, rinse parts thoroughly with water. For spotless drying, rinse with demineralised water.

Respirators (face masks) and protective work equipment have to be cleaned according to the manufacturers specifications.

Composition

TICKOPUR R 33 contains: 5 - 15 % anionic tensides, 5 - 15 % phosphate, <5 % non-ionic tensides, <5 % silicate, complexing agents.

Classification according CLP

Signal word: Danger, GHS05, H315, H318.

regulations concerning sewerage have to be observed.

Completely rinsed jerrycans, bottles and barrels may be recycled via national plastic recycling systems.

Delivery range / dosing aids

		3 3	
1	litre	PE-bottle	Order No. 830
2	litres	PE-handle-bottle	Order No. 883
5	litres	PE-jerrycan	Order No. 831
25	litres	PE-jerrycan	Order No. 835
200	litres	PE-barrel	Order No. 837

Dosing pump for 5- and 25-litres-jerrycan Dosing valves for 25-litres jerrycan

EC material safety data sheet available online as PDF download at: www.dr-h-stamm.info

These specifications are for information only. Subject to technical alterations in composition and change in legal and other safety regulations.

> **DR·H·STAMM** GmbH Chemische Fabrik Heinrichstraße 3 - 4 • 12207 Berlin • GERMANY Tel.: +49 30 76880280 • Fax: +49 30 7734699