

TICKOMED 1

Dental instrument cleaner

for ultrasound • concentrate • mildly alkaline



Cleaning and pre-cleaning of dental instruments and accessories

Concentrate for manual cleaning and pre-cleaning in instrument reprocessing, both in the ultrasonic bath and in the immersion bath process, of

dental, surgical, invasive and non-invasive medical instruments, rotating and non-rotating small and precision instruments as well as small endodontological instruments, accessories and components of medical devices

in dental practices, dental clinics, maxillofacial surgery and orthodontics.

Very high material compatibility, suitable for parts made of metal, and also aluminium and anodised parts, titanium, glass, porcelain, ceramics and plastic. In case of doubt, it must be tested in advance.

Application in the ultrasonic bath (dosage • treatment time)

3 % • 3 – 10 min

Application without ultrasound (dosage • treatment time)

5 % • max. 12 h;

depends on the type and stubbornness of the contamination

5 % • max. 60 min;

for aluminium and its alloys

The application is diluted with water of at least drinking quality and without further additives. Ensure that parts are fully immersed. After treatment, rinse the parts thoroughly with water and dry. Change the working solution immediately if any contamination is visible. Otherwise, replace it at least once each day. Prepare cold for fresh blood and protein residues (protein fixation); solution can be heated up to 60° C for stubborn contamination.

Observe the product information at www.dr-stamm.de.

Only for professional users.

It can be used as a contact liquid (3%) in indirect sonication when using STAMMOPUR Z, STAMMOPUR GR and STAMMOPUR DB in the ultrasonic bath.

Mildly alkaline (pH 9.0 at 1 % in DI water), biodegradable.

Ingredients

5–15% non-ionic tensides, 5–15% cationic tensides, <5% phosphate, complexing agents.

Hazard statements

H315 Causes skin irritation.

H318 Causes serious eye damage.

Precautionary statements

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P305+P351 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if

+P338 present and easy to do. Continue rinsing.

Danger



UFI-Code: NS00-Q0YN-900P-SD06

Serious incidents with the product must be reported to DR.H.STAMM GmbH Chemische Fabrik and to the competent authority.

Basis-UDI-DI: ++G014TM109WP

This instruction for use is a translation. The binding edition is the German version.

S 3090/07e – 2023/01en



DR.H.STAMM GmbH Chemische Fabrik
Heinrichstraße 3-4 • 12207 Berlin (Germany)
Tel.: +49 30768 80 280 • Fax: +49 30 773 46 99
www.dr-stamm.de • info@dr-stamm.de

EC safety data sheets and product information as PDF
download at www.dr-stamm.de