

STAMMOPUR 24

Intensive Cleaning and Disinfection

for Ultrasound • Concentrate

VAH-certified • **MD** C€ 0483

EXAM-Expertise • N-69946

Intensive cleaning and chemical disinfection of medical instruments and accessories

Concentrate for manually disinfecting, non-fixating intensive cleaning and for chemical disinfection, both in the ultrasonic bath and in immersion bath processes for

general, surgical, invasive and non-invasive medical instruments and accessories as well as laboratory equipment, laboratory instruments,

in the CSSD, in the hospital, and in the doctor's practice.

For critical and semi-critical medical devices further reprocessing steps are required; do not use for final disinfection.

Notes on non-medical application

Registered as a biocide under the Biocidal Products Directive (Regulation (EU) No. 528/2012). Also intended for the reprocessing of personal protective equipment (PPE), in particular respiratory masks, protective goggles, and face shields, as well as in the non-medical laboratory for lab equipment and lab instruments.

Very high material compatibility, suitable for parts made of metal, light metal, titanium, glass, porcelain, ceramics, rubber and plastic. In case of doubt, this must be tested in advance.

Application in ultrasonic bath (dosage • treatment time)

1 % • 15 min or 2 % • 5 min:

bactericidal¹, yeasticidal², limited virucidal³ incl. H5N1 and SARS-CoV-2

Application without ultrasound (dosage • treatment time)

1 % • 60 min:

bactericidal¹, yeasticidal²

2 % • 30 min oder 3 % • 15 min:

bactericidal¹, yeasticidal², limited virucidal³ incl. H5N1 and SARS-CoV-2

¹EN 13727, EN 14561, DGHM, high load; ²EN 13624, EN 14562, DGHM, high load; ³EN 14476, high load

The application is carried out diluted with water of at least drinking water quality and without further additives. Ensure that parts are fully immersed. After treatment, rinse the parts thoroughly with water and dry. Change the working solution immediately if there is visible contamination, but not later than after 3 days. Avoid temperatures > 45 °C (protein fixation). Not compatible with aldehyde agents.

Observe the product information at www.dr-stamm.de.

Only for professional users.

Mildly alkaline (pH 9.4 at 1 % in DI water), biodegradable.

100 g contain:

9.9 g bis(3-aminopropyl)dodecylamine, 8.4 g didecylmethylpolyoxyethylammoniumpropionate, < 30 % butane-1,4-diol, 5–15 % butyldiglycol, 5–15 % nonionic tensides, <5 % piperazine, complexing agents, pH regulators.

Hazard statements

- H302 Harmful if swallowed.
H314 Causes severe skin burns and eye damage.
H317 May cause an allergic skin reaction.
H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H336 May cause drowsiness or dizziness.
H361fd Suspected of damaging fertility. Suspected of damaging the unborn child.
H410 Very toxic to aquatic life with long lasting effects.

Precautionary statements

- P280 Wear protective gloves/protective clothing/eye protection/face protection.
P305+P351 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if
+P338 present and easy to do. Continue rinsing.
P308+P313 IF exposed or concerned: Get medical advice/attention.
P310 Immediately call a POISON CENTER/doctor.

Danger



UFI-Code: E800-POU8-F00Q-TAGT

Serious incidents involving the product must be reported to DR.H.STAMM GmbH Chemische Fabrik and the competent authority.

Basis-UDI-DI: ++G014ST2403RG

This instruction for use is a translation. The binding edition is the German version.

S 3460/13e – 2023/02en



DR.H.STAMM GmbH Chemische Fabrik
Heinrichstraße 3-4 • 12207 Berlin (Germany)
Tel.: +49 30768 80 280 • Fax: +49 30 773 46 99
www.dr-stamm.de • info@dr-stamm.de

EC safety data sheets and product information as PDF
download at www.dr-stamm.de