

## Transportwagen und Untergestell für SONOREX Reinigungs- und Spülbäder



Transportwagen SONOREX TECHNIK TW 40  
mit Ultraschallbad SONOREX SUPER RK 1028 CH

### Transportwagen TW - mit feststellbaren Rollen

| Passend für                                                                     | Typ   | Außenmaße          | Arbeitshöhe<br>Oberkante<br>Reinigungswanne | Bestell-<br>Nr. |
|---------------------------------------------------------------------------------|-------|--------------------|---------------------------------------------|-----------------|
| RK 1028 /H, DT 1028 /H<br>RK 1028 CH, DT 1028 CH<br>RM 40.2 /UH /U /H<br>W 65.2 | TW 40 | 550 × 350 × 270 mm | 645 mm<br>745 mm<br>740 mm<br>890 mm        | 8330            |
| RK 1050, DT 1050<br>RK 1050 CH, DT 1050 CH<br>RM 75.2 /UH /U /H                 | TW 75 | 650 × 550 × 250 mm | 650 mm<br>755 mm<br>755 mm                  | 8331            |



### Untergestell UG - mit höhenverstellbaren Füßen, ohne Laufrollen

| Passend für                                                                     | Typ   | Außenmaße              | Arbeitshöhe<br>Oberkante<br>Reinigungswanne                  | Bestell-<br>Nr. |
|---------------------------------------------------------------------------------|-------|------------------------|--------------------------------------------------------------|-----------------|
| RK 1028 /H, DT 1028 /H<br>RK 1028 CH, DT 1028 CH<br>RM 40.2 /UH /U /H<br>W 65.2 | UG 40 | 550 × 350 × 200-220 mm | 575 – 595 mm<br>675 – 695 mm<br>670 – 690 mm<br>820 – 840 mm | 8325            |
| RK 1050, DT 1050<br>RK 1050 CH, DT 1050 CH<br>RM 75.2 /UH /U /H                 | UG 75 | 650 × 550 × 180-200 mm | 580 – 600 mm<br>685 – 705 mm<br>685 – 705 mm                 | 8326            |



## Transport cart and base frame for SONOREX ultrasonic cleaning and rinsing units



transport cart SONOREX TECHNIK TW 40  
with ultrasonic bath SONOREX SUPER RK 1028 CH

### Transport cart TW - with lockable rolls

| For use with                                                                    | Type  | Exterior dimensions | Working height<br>top edge of the<br>cleaning tank | Code-<br>No. |
|---------------------------------------------------------------------------------|-------|---------------------|----------------------------------------------------|--------------|
| RK 1028 /H, DT 1028 /H<br>RK 1028 CH, DT 1028 CH<br>RM 40.2 /UH /U /H<br>W 65.2 | TW 40 | 550 × 350 × 270 mm  | 645 mm<br>745 mm<br>740 mm<br>890 mm               | 8330         |
| RK 1050, DT 1050<br>RK 1050 CH, DT 1050 CH<br>RM 75.2 /UH /U /H                 | TW 75 | 650 × 550 × 250 mm  | 650 mm<br>755 mm<br>755 mm                         | 8331         |



### Base frame UG - with height adjustable feet, without rolls

| For use with                                                                    | Type  | Exterior dimensions    | Working height<br>top edge of the<br>cleaning tank           | Code-<br>No. |
|---------------------------------------------------------------------------------|-------|------------------------|--------------------------------------------------------------|--------------|
| RK 1028 /H, DT 1028 /H<br>RK 1028 CH, DT 1028 CH<br>RM 40.2 /UH /U /H<br>W 65.2 | UG 40 | 550 × 350 × 200-220 mm | 575 – 595 mm<br>675 – 695 mm<br>670 – 690 mm<br>820 – 840 mm | 8325         |
| RK 1050, DT 1050<br>RK 1050 CH, DT 1050 CH<br>RM 75.2 /UH /U /H                 | UG 75 | 650 × 550 × 180-200 mm | 580 – 600 mm<br>685 – 705 mm<br>685 – 705 mm                 | 8326         |



Subject to technical alterations. Dimensions subject to production tolerances.