

## Transportwagen und Untergestell für SONOREX Reinigungs- und Spülbäder



Transportwagen SONOREX TECHNIK TW 40  
mit Ultraschallbad SONOREX SUPER RK 1028 CH

### Transportwagen TW - mit feststellbaren Rollen

Passend für	Typ	Außenmaße	Arbeitshöhe Oberkante Reinigungswanne	Bestell- Nr.
RK 1028 /H, DT 1028 /H RK 1028 CH, DT 1028 CH RM 40.2 /UH /U /H W 65	TW 40	550 × 350 × 270 mm	635 mm 740 mm 740 mm 900 mm	8330
RK 1050, DT 1050 RK 1050 CH, DT 1050 CH RM 75.2 /UH /U /H	TW 75	650 × 550 × 250 mm	645 mm 745 mm 745 mm	8331



### Untergestell UG - mit höhenverstellbaren Füßen, ohne Laufrollen

Passend für	Typ	Außenmaße	Arbeitshöhe Oberkante Reinigungswanne	Bestell- Nr.
RK 1028 /H, DT 1028 /H RK 1028 CH, DT 1028 CH RM 40.2 /UH /U /H W 65	UG 40	550 × 350 × 198-215 mm	560 – 580 mm 665 – 685 mm 665 – 685 mm 825 – 845 mm	8325
RK 1050, DT 1050 RK 1050 CH, DT 1050 CH RM 75.2 /UH /U /H	UG 75	650 × 550 × 175-195 mm	570 – 590 mm 670 – 690 mm 670 – 690 mm	8326



## Transport cart and base frame for SONOREX ultrasonic cleaning and rinsing units



transport cart SONOREX TECHNIK TW 40  
with ultrasonic bath SONOREX SUPER RK 1028 CH

### Transport cart TW - with lockable rolls

For use with	Type	Exterior dimensions	Working height top edge of the cleaning tank	Code- No.
RK 1028 /H, DT 1028 /H RK 1028 CH, DT 1028 CH RM 40.2 /UH /U /H W 65	TW 40	550 × 350 × 270 mm	635 mm 740 mm 740 mm 900 mm	8330
RK 1050, DT 1050 RK 1050 CH, DT 1050 CH RM 75.2 /UH /U /H	TW 75	650 × 550 × 250 mm	645 mm 745 mm 745 mm	8331



### Base frame UG - with height adjustable feet, without rolls

For use with	Type	Exterior dimensions	Working height top edge of the cleaning tank	Code- No.
RK 1028 /H, DT 1028 /H RK 1028 CH, DT 1028 CH RM 40.2 /UH /U /H W 65	UG 40	550 × 350 × 198-215 mm	560 – 580 mm 665 – 685 mm 665 – 685 mm 825 – 845 mm	8325
RK 1050, DT 1050 RK 1050 CH, DT 1050 CH RM 75.2 /UH /U /H	UG 75	650 × 550 × 175-195 mm	570 – 590 mm 670 – 690 mm 670 – 690 mm	8326



Subject to technical alterations. Dimensions subject to production tolerances.